



Multiparameter Indicator for Steam Sterilization



Before exposition.



After exposition.

CHEMGRATOR multiparameter indicators are designed to be used within the packages or specific locations within the sterilization chamber where the sterilizing agent cannot easily access.

The change of color of the multiparameter indicator occurs only in presence of all the critical parameters of a steam sterilization cycle (time, temperature and saturated steam).

Main features of CHEMGRATOR multiparameter indicators:

Hi-quality indicator ink that meets the requirements for ISO 11140-1, EN-867 and IRAM 37101-1 Class 4 standards.

The indicator dot changes from blue to black and the starting and ending colors are printed in each indicator for an immediate interpretation of the results of the sterilization cycle.

The change of color is irreversible. Hence it can be used as a record of the sterilization cycle.

Suitable for regular and flash sterilization cycles.

To be used in steam autoclave at 121°C -7 minutes or 134°C -1 minute.

Shelf Life: 36 months from its date of manufacture